

The Resolution **Revolution**

The Mitaka MM51 is an entirely new microscope design with the highest resolution optics available for super-micro and plastic surgery.



Mitaka MM51

Larger lenses give greater resolution. The larger objective lens of the MM51 combined with Mitaka's advanced 8:1 zoom system creates a microscope with twice the resolution and magnification of other surgical microscopes. This unique design provides a brilliant image at high magnifications, and allows the surgeon to clearly operate on vessels of less than 1 mm. Other companies may add magnification extender accessories; however, they result in a darker image, poor depth of focus, reduced contrast, and degraded image quality.

Mitaka For more information
visit mitakausa.com

改善

A premium microscope
dedicated to addressing the
needs of plastic surgeons.

4K CAMERA AND MONITOR

HIGH RESOLUTION RECORDING DEVICE


ERGONOMIC HAND GRIPS AND FOOT PEDAL

THE ONLY 8:1 ZOOM AVAILABLE FOR SURGERY

FOUR DIFFERENT BINOCULAR HEIGHTS



More than twice the resolution and
the native magnification of any
surgical microscope sold worldwide.



Multispectral, ICG imaging
included with an option to add
fluorescein imaging capability.



**SCHEDULE A
DEMONSTRATION**

Contact Sara Sturgis

sara@mitakausa.com

435.513.0339