



Post-op sinus care should be
simple.



NexPak[®] Intranasal Splint is indicated for use as an intranasal splint intended to minimize bleeding and edema and to prevent adhesions between the septum and the nasal cavity. It is placed in the nasal cavity after surgery or trauma. NexPak is constructed of a patient-comfortable sponge which is manufactured from 100% plant-based polysaccharides.



- Dissolves naturally over 7-14 days or can be easily removed with gentle irrigation and aspiration
- Biocompatible – made from 100% plant-based polysaccharides
- Minimizes bleeding and edema after surgery
- Compressible foam easily conforms to nasal anatomy
- Easy to apply – no mixing or additional prep time necessary
- No special storage conditions required

CAUTION: Federal (U.S.) law restricts this device to sale by or on the order of a physician. NexPak is contraindicated for patients who have known allergies to potato starch. For detailed information regarding indications for use, warnings and precautions, see Instructions for Use.

ORDERING INFORMATION

Part Number:	Description:
9110524	NexPak Intranasal Splint (approx. 0.5" x 1.75")

NexPak is packaged five (5) sterile units per carton. One (1) sterile unit contains one (1) intranasal splint. NexPak is terminally sterilized by gamma irradiation. Store NexPak at room temperature.



To place an order or for more information on NexPak or our other medical technologies, please contact a sales representative by calling 866-612-2568.