



WIGGINS MEDICAL

Home of the Original Black Kerrison

A large, black silhouette of a surgical instrument, specifically a Kerrison rongeur, is positioned in the lower-left quadrant of the page. The instrument is shown from a side profile, highlighting its long, thin handle and the complex, curved headpiece designed for precise bone removal.

Instruments for Spine Surgery

ABOUT WIGGINS MEDICAL

Wiggins Medical has been the leader in specialty surgical instruments since 1992. Our instruments for Spine Surgery are of the highest quality, yet reasonably priced to fit within the hospital's budget.

Wiggins Medical has over 50 years experience in the instrument business. We frequently travel to Germany to meet with the finest German craftsmen to hand-select our instruments and to design new products that are unique to Wiggins Medical.

We were the first company in the United States to offer the black ceramic-coated Kerrisons. Now, there are many copies, but none with the original patented black coating.

In 1995 we introduced the Titanium Distraction Pin into the market.

Compare and see the difference!

After a quick look through the catalog, you will appreciate our offering. Please call our office at (800) 497-0851 for additional information.

Thank you,
Wiggins Medical

TABLE OF CONTENTS

TABLE OF CONTENTS.....	1
WHY BUY OUR INSTRUMENTS	2
PITUITARY RONGEURS	
Basic Patterns - 7", 9", 12", 14"	3
Caspar (fenestrated with teeth) - 7", 9", 14"	4
Curved Foramenotomy - 7"	5
Decker - 6", 10"	5
Pituitary with tooth - 7"	6
Peapods - 7"	6
KERRISON RONGEURS	
Spurling (large handle) (16mm opening) - 8"	7
Kerrisons - 7", 8", 9", 10", 12", 14"	8
Kerrisons - 7", 8", 9", 10", 12", 14"	9
Down Biters - 8"	10
Curved Foramenotomy - 8"	11
Bayonet - 7" (with 7" working length)	12
Rotating Shaft - 8", 9", bayonet.....	13
Endoscopic rotating shaft - 10", 12", 14"	14
Detachable - 8", 9"	15
TITANIUM DISTRACTION PINS.....	16
CONTACT INFORMATION AND TERMS OF SALE	17

All Instruments have the Patented Black Ceramic Coating

WHY BUY OUR KERRISONS ?

We offer a wide variety of kerrison patterns and styles. We have them in a standard handle, large spurling handle, with and without ejector, standard, thin and ultra-thin footplates, straight, bayonet, curved, detachable, rotating shafts, and many other new designs. All have our PATENTED BLACK CERAMIC COATING. If you do not see the length or pattern you need, please contact us at (800) 497-0851 for a quotation.

Wiggins Medical introduced the ceramic-coated design to the US market in July 1995. Our intent was to provide a bone-biting instrument with a razor-smooth cutting action, produce less hand fatigue for the surgeon, and reduce the problem of glare.

An unexpected result was that the patented special coating actually extended the cutting life of the instrument. We have some customers that have never had to have their instruments resharpended.

We can say, with confidence, that we have the "Thinnest Thin Footplate", even in a 1mm size. Our thin footplate pattern should only be used in cervical or other delicate procedures.

If screws have ever fallen out of your kerrisons...our kerrison screws are reverse-threaded so they will not back out. As you squeeze the handle, the screws stay securely in place.

WHY BUY OUR PITUITARY RONGEURS ?

Standard Pituitary Rongeurs - We were the first to offer an ultra-thin, 1mm, micro intervertebral disc rongeur. Of course, the smoothness, anti-glare, and precise cutting edge, are a result of fine German craftsmanship and our black ceramic coating.

Fenestrated Pituitary Rongeurs - This is a very popular style. In our version, the cup shape is slightly tapered to allow better access to the disc space. There are teeth placed proximally in the cup to hold tissue better. The fenestrations allow tissue to expand into the cup and give the surgeon a more precise bite.

WHY BUY OUR TITANIUM DISTRACTION PINS VS. STAINLESS STEEL ?

- | | |
|-------------------------------------|---|
| 1. Lighter | 6. Self Tapping Tip (Most pre-drilling not needed) |
| 2. Stronger | 7. Shoulder fits Close to the Bone (Greater Lateral Stability and More Secure Hold) |
| 3. Safer | 8. Greater Imaging Quality |
| 4. Implant Grade | 9. Pre-packaged Sterile, Single Use |
| 5. Nickel-Free (No Nickel Allergen) | 10. We Never backorder |

PITUITARY RONGEURS



STRAIGHT

	7"	9"	14"
.5mm	52-051-05B		
1mm	52-051-10B	52-110-10B	
1.5mm	52-051-15B	52-110-15B	
2mm	52-051-20B	52-110-20B	52-119-20/36B
3mm	52-051-30B	52-110-30B	52-119-30/36B
4mm	52-051-40B	52-110-40B	52-119-40/36B
5mm	52-051-50B	52-110-50B	
6mm		52-110-60B	



UP

	7"	9"	14"
.5mm	52-053-05B		
1mm	52-053-10B	52-111-10B	
1.5mm	52-053-15B	52-111-15B	
2mm	52-053-20B	52-111-20B	52-120-20/36B
3mm	52-053-30B	52-111-30B	52-120-30/36B
4mm	52-053-40B	52-111-40B	52-120-40/36B
5mm	52-053-50B	52-111-50B	
6mm		52-111-60B	



DOWN

	7"	9"	14"
2mm	52-055-20B	52-046-20/23B	52-046-20/36B
3mm	52-055-30B	52-046-30/23B	52-046-30/36B
4mm	52-055-40B	52-046-40/23B	52-046-40/36B
5mm	52-055-50B		



12" WITH TEETH

	Straight	Up
2mm	52-051-20/30SB	52-053-20/30SB
3mm	52-051-30/30SB	52-053-30/30SB
4mm	52-051-40/30SB	52-053-40/30SB

CASPAR FENESTRATED (with teeth) Pituitary Rongeurs



STRAIGHT

	7"	9"	14"
2mm	52-096-20B	52-096-20/23B	52-096-20/36B
3mm	52-096-30B	52-096-30/23B	52-096-30/36B
4mm	52-096-40B	52-096-40/23B	52-096-40/36B
5mm	52-096-50B	52-096-50/23B	52-096-50/36B
6mm	52-096-60B	52-096-60/23B	



UP

	7"	9"	14"
2mm	52-097-20B	52-097-20/23B	52-097-20/36B
3mm	52-097-30B	52-097-30/23B	52-097-30/36B
4mm	52-097-40B	52-097-40/23B	52-097-40/36B
5mm	52-097-50B	52-097-50/23B	
6mm	52-097-60B	52-097-60/23B	




DOWN

	7"	9"	14"
2mm	52-106-20B	52-106-20/23B	52-106-20/36B
3mm	52-106-30B	52-106-30/23B	52-106-30/36B
4mm	52-106-40B	52-106-40/23B	52-106-40/36B
5mm	52-106-50B	52-106-50/23B	
6mm	52-106-60B	52-106-60/23B	

CURVED (FORAMENOTOMY) PITUITARY

7"
 1mm 52-008-10/A/B
 2mm 52-008-20/A/B
 3mm 52-008-30/A/B



DECKER 2mm X 6mm CUP

6"
 2mm 51-011-15B

10"
 1mm 51-011-10/25B (1mm X 6mm CUP)
 2mm 51-011-20/25B (2mm X 6mm CUP)
 3mm 51-011-30/25B (3mm X 6mm CUP)



PEAPOD

Up Biting

7"

1mm	52-068-10B
1.5mm	52-068-15B
2mm	52-068-20B
2.5mm	52-068-25B
3mm	52-068-30B
4mm	52-068-40B



STRAIGHT PEAPOD

Straight

7"

1mm	52-068-10/GB
1.5mm	52-068-15/GB
2mm	52-068-20/GB
3mm	52-068-30/GB



PITUITARY WITH TOOTH

Pituitary with Tooth

7"

2mm	52-000-02/19B
3mm	52-000-03/19B



SPURLING (Large Handle) KERRISON 40° Up, Standard Footplate



8" (40° up) no Ejector

2mm	52-262-02B
3mm	52-262-03B
4mm	52-262-04B
5mm	52-262-05B
6mm	52-262-06B

8" (40° up) with Ejector

52-262-02EB
52-262-03EB
52-262-04EB
52-262-05EB
52-262-06EB

8" (90° up) no Ejector

2mm	52-242-02B
3mm	52-242-03B
4mm	52-242-04B

8" (40° up) with Ejector Detachable Version

2mm	52-762-02EB
3mm	52-762-03EB
4mm	52-762-04EB
5mm	52-762-05EB

KERRISONS

40° Up, 9mm Opening



Standard Footplate

Thin Footplate

Ultra-Thin Footplate

Standard Footplate with Ejector

	7"	8"	9"	10"	12"	14"
2mm	52-161-02EB	52-162-02EB	52-165-02EB	52-166-02EB	52-168-02EB	52-168-02/36EB
2.5mm		52-162-025EB				
3mm	52-161-03EB	52-162-03EB	52-165-03EB	52-166-03EB	52-168-03EB	52-168-03/36EB
4mm	52-161-04EB	52-162-04EB	52-165-04EB	52-166-04EB	52-168-04EB	52-168-04/36EB
5mm	52-161-05EB	52-162-05EB	52-165-05EB	52-166-05EB	52-168-05EB	52-168-05/36EB
6mm		52-162-06EB	52-165-06EB	52-166-06EB	52-168-06EB	52-168-06/36EB

Standard Footplate no Ejector

	7"	8"	9"	10"	12"	14"
.6mm	52-161-006B	52-162-006B				
1mm	52-161-01B	52-162-01B	52-165-01B	52-166-01B	52-168-01B	
1.5mm		52-162-015B				
2mm	52-161-02B	52-162-02B	52-165-02B	52-166-02B	52-168-02B	52-168-02/36B
2.5mm		52-162-025B				
3mm	52-161-03B	52-162-03B	52-165-03B	52-166-03B	52-168-03B	52-168-03/36B
4mm	52-161-04B	52-162-04B	52-165-04B	52-166-04B	52-168-04B	52-168-04/36B
5mm	52-161-05B	52-162-05B	52-165-05B	52-166-05B	52-168-05B	
6mm		52-162-06B	52-165-06B		52-168-06B	

Thin Footplate with Ejector

	7"	8"	9"	10"	12"
2mm	52-181-02EB	52-182-02EB	52-185-02EB	52-186-02EB	52-188-02EB
2.5mm		52-182-025EB			
3mm	52-181-03EB	52-182-03EB	52-185-03EB	52-186-03EB	52-188-03EB
4mm	52-181-04EB	52-182-04EB	52-185-04EB	52-186-04EB	52-188-04EB
5mm	52-181-05EB	52-182-05EB	52-185-05EB	52-186-05EB	52-188-05EB
6mm		52-182-06EB	52-185-06EB		

Thin Footplate no Ejector

	7"	8"	9"	10"	12"
1mm	52-181-01B	52-182-01B	52-185-01B	52-186-01B	52-188-01B
1.5mm		52-182-015B			
2mm	52-181-02B	52-182-02B	52-185-02B	52-186-02B	52-188-02B
3mm	52-181-03B	52-182-03B	52-185-03B	52-186-03B	52-188-03B
4mm	52-181-04B	52-182-04B	52-185-04B	52-186-04B	52-188-04B
5mm	52-181-05B	52-182-05B	52-185-05B	52-186-05B	52-188-05B
6mm		52-182-06B	52-185-06B		

Ultra Thin Footplate with Ejector

	8"
2mm	52-212-02EB
3mm	52-212-03EB
4mm	52-212-04EB


Ultra Thin Footplate with No Ejector

	8"
1mm	52-212-01B
2mm	52-212-02B
3mm	52-212-03B
4mm	52-212-04B

FOOTPLATE COMPARISON CHART

Type	Jaw Width	Thickness (*)
Standard Footplate Kerrison	1mm	2.1mm +/- 0.1mm
	2mm	2.3mm +/- 0.1mm
	2.5mm	2.5mm +/- 0.1mm
	3mm	2.6mm +/- 0.1mm
	3.5mm	2.7mm +/- 0.1mm
	4mm	2.8mm +/- 0.1mm
	5mm	3.0mm +/- 0.1mm
Thin Footplate Kerrison	6mm	3.3mm +/- 0.1mm
	1mm	1.7mm +/- 0.1mm
	2mm	1.8mm +/- 0.1mm
	2.5mm	1.9mm +/- 0.1mm
	3mm	2.0mm +/- 0.1mm
	4mm	2.2mm +/- 0.1mm
Ultra-Thin Footplate Kerrison	5mm	2.4mm +/- 0.1mm
	1mm	1.3mm +/- 0.1mm
	2mm	1.5mm +/- 0.1mm
	3mm	1.6mm +/- 0.1mm
	4mm	1.8mm +/- 0.1mm
	5mm	2.0mm +/- 0.1mm

(*) Thickness is measured at the lower end of cut at 90 degree angle towards cut



DOWN BITERS (40° + 90°)



90° down Standard Footplate

8"
no Ejector

2mm	52-152-02B
3mm	52-152-03B
4mm	52-152-04B

40° down Standard Footplate

8"
no Ejector

2mm	52-172-02B
3mm	52-172-03B
4mm	52-172-04B

90° down Thin Footplate

8"
no Ejector

1mm	52-152-01PL/B
2mm	52-152-02PL/B
3mm	52-152-03PL/B

40° down ThinFootplate

8"
no Ejector

1mm	52-172-01PL/B
2mm	52-172-02PL/B
3mm	52-172-03PL/B

CURVED FORAMENOTOMY KERRISON



Thin Footplate


8"		8"		8"	
Standard Curve		Slight Curve w/ Short Footplate		Slight Curve	
2mm	52-300-02B	2mm	52-305-02B	2mm	52-350-02B
3mm	52-300-03B	3mm	52-305-03B	3mm	52-350-03B
4mm	52-300-04B	4mm	52-305-04B	4mm	52-350-04B

Ultra-Thin Footplate

8"		8"		8"	
Standard Curve		Slight Curve w/ Short Footplate		Slight Curve	
2mm	52-300-02SPFB	2mm	52-305-02SPFB	2mm	52-350-02SPFB
3mm	52-300-03SPFB	3mm	52-305-03SPFB	3mm	52-350-03SPFB
4mm	52-300-04SPFB	4mm	52-305-04SPFB	4mm	52-350-04SPFB

BAYONET KERRISON

40° Up, Straight



Standard Footplate

7" no Ejector		7" with Ejector	
1mm	52-610-01B	2mm	52-610-02EB
2mm	52-610-02B	3mm	52-610-03EB
3mm	52-610-03B	4mm	52-610-04EB
4mm	52-610-04B	5mm	52-610-05EB

Thin Footplate

7" no Ejector		7" with Ejector	
1mm	52-616-01B	2mm	52-616-02EB
2mm	52-616-02B	3mm	52-616-03EB
3mm	52-616-03B	4mm	52-616-04EB
4mm	52-616-04B	5mm	52-616-05EB

ROTATING SHAFT KERRISON

40° Up, 14mm opening

14 mm Opening



Shafts with Thin Footplates

8" Shaft with Ejector

1mm	52-482-01B
2mm	52-482-02EB
3mm	52-482-03EB
4mm	52-482-04EB
5mm	52-482-05EB
6mm	52-482-06EB

9" Shaft with Ejector

1mm	52-485-01B
2mm	52-485-02EB
3mm	52-485-03EB
4mm	52-485-04EB
5mm	52-485-05EB

7" Bayonet Shaft with Ejector

1mm	52-416-01B(no ejector in 1mm)
2mm	52-416-02EB
3mm	52-416-03EB
4mm	52-416-04EB

Handle Part Number 52-400-00B



9" STRAIGHT SHAFT



8" STRAIGHT SHAFT




7" WORKING DISTANCE
BAYONET SHAFT



HANDLE IS FOR USE WITH ALL ROTATING SHAFTS
52-400-00B

ROTATING ENDOSCOPIC KERRISON SHAFT

40° Up, No Ejector



3mm tube and 3mm tip diameter 52-439-03B



3.5mm tube and 3.5mm tip diameter 52-439-035B



4mm tube and 4mm tip diameter 52-439-04B



52-439-04B



52-439-035B



52-439-03B



**HANDLE
FOR USE WITH
ALL SHAFTS**
52-400-00B

Available in Thin Footplate only

	10"		12"		14"
3mm	52-437-03B	3mm	52-438-03B	3mm	52-439-03B
3.5mm	52-437-035B	3.5mm	52-438-035B	3.5mm	52-439-035B
4mm	52-437-04B	4mm	52-438-04B	4mm	52-439-04B

“DETACHABLE” KERRISON 40° Up with Ejector



14mm opening



Curved Detachable Kerrison Thin Footplate



8"

Standard Curved Tip with Ejector 9mm Opening

- 2mm 52-742-02ECB
- 3mm 52-742-03ECB

Thin Footplate with Ejector

8"

9"

- | | | | |
|-------|-------------|-------------|-------------------|
| 1mm | 52-742-01B | 52-745-01B | No Ejector on 1mm |
| 1.5mm | 52-742-015B | | |
| 2mm | 52-742-02EB | 52-745-02EB | |
| 3mm | 52-742-03EB | 52-745-03EB | |
| 4mm | 52-742-04EB | 52-745-04EB | |
| 5mm | 52-742-05EB | 52-745-05EB | |
| 6mm | 52-742-06EB | 52-745-06EB | |

Standard Footplate with Ejector

8"

- 1mm 52-732-01B
- 2mm 52-732-02EB
- 3mm 52-732-03EB
- 4mm 52-732-04EB
- 5mm 52-732-05EB

Large Handle Standard Footplate with Ejector

8"

- 2mm 52-762-02EB
- 3mm 52-762-03EB
- 4mm 52-762-04EB
- 5mm 52-762-05EB
- 6mm 52-762-06EB



8"

Slightly Curved Tip with Ejector 9mm Opening

- 2mm 52-742-02ELCB
- 3mm 52-742-03ELCB

Avoid Nickel Allergy

DISTRACTION PINS TITANIUM VS STEEL

WIGGINS



Titanium

**Flush Mount Collar for
greater lateral stability**

No Halo with Imaging

**No Nickel for Potential
Allergic Reaction**

**Self Tapping Eliminates
most needs for pre-drilling**

Implant Grade

WE NEVER BACKORDER

THEIRS



Steel

NO

Visible Halo with Imaging

**All Steel has Nickel for
Potential Allergic Reaction**

NO

NO

NO

Size	Item#	Quantity/Box
12mm	67-012	5/Box
12mm	67-012-10	10/Box
14mm	67-014	5/Box
14mm	67-014-10	10/Box
16mm	67-016	5/Box
16mm	67-016-10	10/Box

CONTACT INFORMATION AND TERMS OF SALE

PURCHASE ORDERS MAY BE SENT DIRECTLY TO:

WIGGINS MEDICAL
818 Ellison Avenue
Cincinnati, Ohio 45226
www.wigginsmedical.com

Fax: (850) 942-6595
Customer Service: (800) 497-0851
info@wigginsmedical.com

PRICES

Prices are listed in Wiggins Medical current price list and are subject to change without notice. Freight is FOB destination, prepaid and added to the invoice.

WARRANTY

All instruments are warranted by the manufacturer against defects in workmanship and materials. If the products are not used in accordance with the intended purpose, altered, repaired or modified by persons not authorized by Wiggins Medical, the warranty will be void. Wiggins Medical's sole obligation and customer's exclusive remedy for breach of any warranty, shall be at Wiggins Medical's option, to repair or replace the product. There are no other expressed or implied warranties, including any warranty of merchantability or fitness for a particular purpose.

REPAIR PROGRAM

Wiggins Medical offers a full-service, warranted repair program. Check with customer service for details.

CUSTOM INSTRUMENTS

Wiggins Medical will discuss the feasibility of a custom instrument. We are happy to quote any customization that is deemed possible to manufacture.

RETURN GOODS POLICY

1. All returns must be authorized by Wiggins Medical.
2. Returns must be in resalable condition within 30 days of the invoice date. A 15% reprocessing fee will be charged unless the return is due to Wiggins Medical error.
3. All returns are for credit only.
4. All returns must have an RGA number issued by the Wiggins Medical customer service department.

NON-RESALABLE MERCHANDISE INCLUDES THE FOLLOWING:

Products not purchased from Wiggins Medical, products that have been sterilized and/or used, products attempted to return 30 days beyond the invoice date, sterile products opened or removed from their original container, custom made or discounted products, and products not in the original packaging.



Visit **Wiggins Medical** on the web at
www.wigginsmedical.com

818 Ellison Avenue
Cincinnati, Ohio 45226

Phone: (800) 497-0851
Fax: (850) 942-6595
info@wigginsmedical.com